



PREVOR

ANTICIPATE AND SAVE

Toxicology Laboratory & Chemical Risk Management

**Choose
Comfort**

DRESSILK[®]

Temporary skin
So soft on you

For superficial wounds
Burns or skin graft donor sites

**Excellent compromise
Cost/effectiveness**



Perfect for superficial wounds
Burns or skin graft donor sites



**100%
natural silk**

1 Simple & easy to use

In just 1 application Dressilk®
barrier protects until wound
closure (30 days maximum)



2 Comfortable & user-friendly

Soft, it facilitates application
due to its enlargement potential
in 2D



MORE INFORMATION ON WWW.DRESSILK.COM

PREVOR & DRESSILK® - Trademarks, names and logo types registered by PREVOR. © 01/2018 PREVOR

3 Excellent compromise

Dressilk® provides the best environment for wound healing while limiting both patient's pain (no or less dressing changes) and hospital care costs

DRESSILK® is removed easily
from any part of the body
after re-epithelialization

36 YEARS OLD PATIENT

1ST DAY

Wound before and after debridement




Application of Dressilk®



Application of secondary dressing



8TH DAY Secondary dressing removed, wound is dry



15TH DAY Re-epithelialization continues, Dressilk® may be trimmed slowly



29TH DAY Complete wound healing



Why choose DRESSILK®?

For staff

- **Simple and user-friendly**
- **Easy** to apply
- **Remains** in place until complete healing
- **Removed** easily
- **Good healing** at short term and **good scarring** at long term

For patients

- **Comfortable** in-use
- **Remains** in place until complete healing
- **Reduces pain** as no or reduced dressing changes needed
- **Good scarring**

For hospitals

- **Saves personnel's time**
- **Good healing**
- **Reduces costs**

Technical information

Medical Device, Class IIb
CE 0459
ARTG n° 281295

- 100 % natural silk
- Enlargement potential in 2 dimensions
- Available in 2 sizes: 10 x 10 cm & 20 x 30 cm
- Allows wound observation in the first days
- Perfect for any superficial wounds, burns or skin graft donor sites
- Satisfactory healing time
- Limits inflammation and infection
- Does not disintegrate
- Excellent compromise between cost and efficiency
- Non irritant, non cytotoxic and non sensitizing



Ref: Schulz A, Depner C, Lefering R, Kricheldorf J, Kästner S, Fuchs PC, Demir E. A prospective clinical trial comparing Biobrane® Dressilk® and PolyMem® dressings on partial-thickness skin graft donor sites. Burns. 2015 Feb 23. pii: S0305-4179(14)00450-1. doi: 10.1016/j.burns.2014.12.016

Ref: Schiefer JL, Arens E, Grigutsch D, Rath R, Hoffmann A, Fuchs PC, Schulz A. A prospective intra-individual evaluation of silk compared to Biobrane for the treatment of superficial burns of the hand and face. Burns. 2016 Dec 30. pii: S0305-4179(16)30354-0. doi:10.1016/j.burns.2016.09.005.

Ref: Schiefer JL, Rath R, Ahrens E, Grigutsch D, Gräff I, Stromps JP, Fuchs PC, Schulz A. Does Dressilk® represent an alternative to Biobrane®? Evaluation of scar quality after treatment of superficial burns of hands and faces - an intra-individual comparison. Burns. 2017 Aug 30. pii: S0305-4179(17)30418-7. doi: 10.1016/j.burns.2017.07.024.

DRESSILK®_2018_01_Flyer_GB_1

Design : Planète Graphique Studio - 01 42 67 67 80 - Paris 17



PREVOR

ANTICIPATE AND SAVE

Toxicology Laboratory & Chemical Risk Management

Moulin de Verville - 95760 VALMONDOIS - France
Tél. : +33 (0)1 30 34 76 76 - Fax : + 33 (0)1 30 34 76 70
mail@prevor.com - www.prevor.com